



High Support Car Seat

Start 2 Model





NECK CUSHION

Fits around the neck to support the head.



SOUND SYSTEM

Listen to music through the integral speakers in the headrest, comes with practical MP3 player pouch on the seat.



FLEXIBLE TRUNK SUPPORTS

Securely fastening around the child to provide additional lateral stability. A waistcoat harness is also available.



TURNTABLE

Allows easier entrance and exit from vehicles, using a straight forward mechanism with smooth action.



We have added extra supportive features tailoring the chair to the requirements of children with special needs.

Adjustments in height, width and seat length allows the seat to grow with the child whilst providing an ergonomic and comfortable seat position. The wide side cushions can be adjusted individually and independently from the height setting.

Security comes in the form of a stable, aluminium frame construction which absorbs energy in a crash, whilst a specially formed ramp prevents the child from slipping under the safety belt.

The seat has passed crash testing to standard ECE-R44/04 (2-3 grouping).



DIMENSIONS (mm)	MIN	MAX
Backrest Height	560	760
Seat Width	-	330
Seat Depth	300	360
User age / Weight range Approx 3-12 Years / 15-36kg		

Tilt adjustment enables a recline position (only for weight class II, approx 3-6 years, 15 – 25 kg).

All covers are breathable, removable and machine washable at 30°C. The seat is also foldable for transport.

www.jcmseating.co.uk

JCM Seating Solutions
15-18 Maxwell Road, Woodston Industrial Estate
Peterborough Cambridgeshire PE2 7HU

T: 01733 405830 E: enquiries@jcmseating.co.uk
F: 01733 405838 W: www.jcmseating.co.uk

