



High Support Car Seat

Young Sport Model





CARRY HANDLE



5-POINT HARNESS



INTEGRATED TRUNK SUPPORT



ADJUSTABLE HEADREST



TURNTABLE



The car seat has passed crash testing to standard ECE-R44/04 (1-3 grouping)

Extremely long life span - Suitable for children from approx 9 months -12 years within a weight range of 9-36kg.

The seat comes with integrated trunk support as well as an enhanced trunk support option if required. All covers are removable and can be machine-washed at 30°C

The 5-point harness is for a child of group 1 (9-18kg) and can be adjusted to fit.

The seat can be easily fitted by setting into the car and secured with the cars own 3-point-belt.

Our turntable mechanism option also allows effortless entrance and exit to vehicles.

DIMENSIONS (mm)	MIN	MAX
Backrest Height	530	640
Seat Width	185	285
Seat Depth	300	
User Age / Weight Range: Approx 9 months-12 years / 9-36kg		



NECK ROLL OPTION



ENHANCED TRUNK SUPPORT OPTION

www.jcmseating.co.uk

JCM Seating Solutions
15-18 Maxwell Road, Woodston Industrial Estate
Peterborough Cambridgeshire PE2 7HU

T: 01733 405830 E: enquiries@jcmseating.co.uk
F: 01733 405838 W: www.jcmseating.co.uk

