



Vision

standing frame
with power lift





LOWERED HANDRAILS

Ensures effortless wheelchair to upright standing.



PADDED KNEE PADS

One-handed height adjustment of the knee pads.



POSTERIOR SUPPORT

Can be folded down to allow easier entrance and exit.



REMOTE CONTROL

Comes with a cordless, easy to use remote control.

The Vision standing frame with power lift offers optimum aid to meet the complete requirements of the user.

The Vision is suitable for conditions such as cerebral palsy, paraplegia, tetraplegia, multiple sclerosis and many others.

Close to floor entry makes standing from a wheelchair straightforward with minimal effort required from the user.


Four pivoted wheels with brakes which lock fast enhance stability and allow for easy manoeuvrability.



Features

- Pressure relieving posterior harness.
- Depth & tilt alterable activity table with user friendly, spring-assisted height adjustment.
- Posterior support which can be folded down.
- Depth, width, height and angle adjustable padded knee supports.
- Depth and width adjustable heel rests.
- Four pivoted wheels with lock-fast brakes.

BASIC MODEL - DIMENSIONS (cm)	
Foot Board - Table	93 - 120
Foot Board - Knee Supports	38 - 56
Height of Handrails (from foot board)	80 - 104
Floor - Upper edge of foot board	8
Max. Length	110
Base Width	78
Table (width x depth)	75 x 60
Exercise Tension	230 v
Weight	62 kg
Max. Load	150 kg



Accessories

- Enlarged belt system
- Posterior Support
- Back support with headrest

www.jcmseating.co.uk

JCM Seating Solutions
15-18 Maxwell Road, Woodston Industrial Estate
Peterborough Cambridgeshire PE2 7HU

T: 01733 405830 **E:** enquiries@jcmseating.co.uk
F: 01733 405838 **W:** www.jcmseating.co.uk

